



PM40 Patient Monitor

Making Parameters Monitoring Simple

Simple clinical operation procedures



PM40 series have customized parameter layout and big numeric display, it's much closer to clinical application.



It has various hotkey editing function and convenient for faster operation. The single-interface with three stage menu design ensures 95% of operations can be completed within 2 steps.



The multiple display interfaces, including normal, big numeric, full screen 7-lead.

Excellent Performance

- ✓ Compact design with hidden handle, easy to carry
- ✓ Standard wired network to connect CMS
- ✓ Intelligent blood pressure measurement technology
- ✓ Optional highly sensitive touch screen offering smooth and easy operation
- ✓ 12.1-inch screen, high resolution 1280*800
- ✓ Standard visual/audio alarm function
- ✓ Standard inbuilt lithium battery for 2 hrs
- ✓ Dynamic EWS score, drug concentration calculation & Hemodynamic monitoring



Accurate Measurement



SpO2

- Anti-interference SpO2 algorithm provides accurate measurement
- Excellent weak perfusion (reach 0.075%)

NIBP

- Dual-path pressure sensor, dual MCUs, more secure and reliable
- Intelligent algorithm patented has the advantages of short inflation time, fast measurement and less patient resistance

ECG

- High T-wave suppression capability (2.4 mV), solve the issue of doubled HR
- High isolation voltage, 4000 Vrms, providing greater safety and reliability

IBP

- IBP waveform superposition function facilitates the IBP comparison and analysis of different parts
- Large Range: -50 - 400mmHg

RESP

- More respiratory lead select: RA-LA、RA-LL、LA-RL、LL-RL

CO2

- CO2 module Low flow rate, minimum flow rate 50ml/min
- Airway temperature and humidity measurement and compensation
- Ultra-long life gas sampling tube, up to 10 days or more
- Ultra-long life brush motor (More than 30,000 hours)

VERITAS MEDTECH PRIVATE LIMITED

S.No. 108, Sri Vignesh Nagar, Kumudham Nagar, Cheran ma Nagar, Villankurichi, Coimbatore, Tamil Nadu 641035



+91 93440 32972
+91 95665 32972



contactus@veritasmtech.com



www.veritasmtech.com



MD-42 : TN/CBN/MD42/00044



GSTIN : 33AALCV4506EIZI