



DEXUP Patient Warmer

Maintains perioperative normothermia to minimize hypothermia-related risks. Integrated safety sensors and a lightweight, trolley-mounted design ensure efficient and reliable patient warming.

◆ Certifications



◆ Specifications

Parameter	Specification
Model	DEXUP PW 400
Dimensions (mm)	250 (L) * 300 (W) * 340 (H) (not including extrusions)
Weight	~5.5 kg
Air Flow	Accuracy: ±10 m3/h Level 0: 15 m3/h Level 1: 30 m3/h Level 2: 45 m3/h Level 3: 60 m3/h Level 4: 70 m3/h Level 5: 85 m3/h
Maximum Output Pressure	>2600Pa
Temperature Range	33°C to 43°C at 0.5°C increments; accuracy ±1.5°C
Modes of Operation	Ambient Mode, Warming Mode
Alarms	Overheat alarm, Temperature control alarm, High temperature alarm, Low temperature alarm, Heater error alarm, Fan error alarm, Temperature sensor alarm
Power Supply	a.c.220V, 50Hz, Power: 1800VA(Max.)
Safety Class	Class I Type BF, IPX1 (Anti dripping sealing unit)
Temperature	+10°C to +40°C
Humidity	30% to 75%
Atmospheric Pressure	700hPa to 1060hPa
Filter	HEPA Filter 0.3 micron / Foam Filter

◆ Features



**Precise
Temperature Control**



**High Speed Heat
Circulation**



**Versatile Blanket
Compatibility**



**Easy Cleaning &
Maintenance**



**Lightweight
Portable Design**



**Advanced
Safety Features**



**User Friendly
Control Panel**



**Durable & High
Performance Build**



◆ Areas of Application

- 01 Pre And Post Surgical Areas
- 02 Labour Rooms
- 03 Intensive Care Units
- 04 High Dependency Units
- 05 Emergency And Trauma Care

◆ Advantages



Hepa Air Filtration

Ensures clean air supply



Built In Overheat Sensors

Detect overheating with Alarm indicator.



Compact Design

Lightweight, portable design



Efficient Heat Transfer

Optimum heat transfer and delivery towards patient.



Trolley Compatible


Convenient trolley mounting option.




Undisturbed Airflow

No disruption of airflow in the surrounding environment.


◆ Blankets



Full Body Blanket



Lower Body Blanket



Upper Body Blanket

VERITAS MEDTECH PRIVATE LIMITED

S.No. 108, Sri Vignesh Nagar, Kumudham Nagar, Cheran ma Nagar, Villankurichi, Coimbatore, Tamil Nadu 641035



+91 93440 32972
+91 95665 32972



contactus@veritasmtech.com



www.veritasmtech.com



MD-42 : TN/CBN/MD42/00044



GSTIN : 33AALCV4506EIZI