



CE

SYRINGE PUMP



INFUSION PUMP

ISO
13485:2016

CE



SYRINGE PUMP Model : SIP-201



Model : SIP-201

- ✓ Accepts all standard syringes 5ml, 10ml, 20ml/30ml, 50ml/60ml.
- ✓ Large COLOR TFT display with adjustable brightness.
- ✓ Five work modes: Rate Mode, Time Mode, Dose Mode, Sequential Mode and Drug Library Mode (Optional).
- ✓ Manual and Automatic Bolus function. KVO and Bolus functions.
- ✓ Anti-bolus system to reduce significantly after occlusion
- ✓ Identifies the specifications of syringe automatically.
- ✓ 6 Hours Internal Battery continuous operation.
- ✓ Audible and visual alarm: Syringe Holder Loose, Syringe Plunger Disengagement, Syringe Not Calibrated, Occlusion, Syringe Empty, Infusion Complete, Syringe Near Empty, Battery Low
- ✓ Accuracy $\pm 2\%$ with correct calibration
- ✓ Total Volume Range 0.1-9999.9ml (0.1 ml step) 0.1ml/h-2000.0ml/h (0.1ml/h or 1ml/h step) TIME 00:01-99:59min (1min step).
- ✓ Bolus Rate 5ml: 0.1-200ml/h, 10ml:0.1-400ml/h, 20ml:0.1-800ml/h, 30ml:0.1-1200ml/h, 50ml:0.1-2000ml/h

INFUSION PUMP Model : IVP-302

- ✓ Compatible with all standard IV sets.
- ✓ Large COLOR TFT display with adjustable brightness.
- ✓ Five work modes: Rate Mode, Time Mode, Dose Mode, Sequential Mode and Drug Library Mode. Drug library with 1030 drug names (Optional).
- ✓ Manual/Automatic Bolus. & Anti-bolus system. KVO and Bolus functions.
- ✓ 6 Hours Internal Battery continuous operation.
- ✓ Visual and audible Alarm system for all failures and Hi - Low flow
- ✓ Dynamically display working status.
- ✓ Three levels of the occlusion: high, middle, and low.
- ✓ Accuracy $\pm 2\%$
- ✓ Delivery rate range 0.1ml/h 1200ml/h (0.1ml/h 99.9ml/h, 0.1ml/h step; >99.9ml/h, 1ml/h step)
- ✓ Delivery time range 00h01min 99h59min (1 min step)
- ✓ Volume limit range: 0.0-9999ml (1 ml step)
- ✓ Total Volume Range 0.0-9999ml (1 ml step)
- ✓ Purge Rate 600.0ml/h (Adjustable) KVO Rate 0.1-5.0ml/h. Programmable and adjustable



Model : IVP-302

VERITAS MEDTECH PRIVATE LIMITED

S.No. 108, Sri Vignesh Nagar, Kumudham Nagar, Cheran ma Nagar, Villankurichi, Coimbatore, Tamil Nadu 641035



+91 93440 32972
+91 95665 32972



contactus@veritasmtech.com



www.veritasmtech.com



MD-42 : TN/CBN/MD42/00044



GSTIN : 33AALCV4506EIZI